


FROM MOLEX SHANNON
DATE SENT : _____

FROM MOLEX SHANNON
DATE SENT : _____

FROM MOLEX SHANNON
DATE SENT : _____

FROM MOLEX SHANNON
DATE SENT : _____


FROM MOLEX SHANNON
DATE SENT : _____



BILL OF MATERIALS

| NO. | QUANTITY | DESCRIPTION |
|-----|----------|---|
| ▽ | 1 | 50 RG 178 / 196 OHM COAXIAL CABLE PER MIL. C. IT STANDARD |
| ▽2 | 1 | STRAIGHT BNC PLUG CONNECTOR FOR RG178/196 CABLE |
| ▽3 | 1 | RIGHT ANGLED RF PLUG PER SDA-91205 |

WIRING DIAGRAM : COAXIAL CABLE ASSEMBLY



NOTES:

1. PACKAGING: LOOSE IN BAGS.

ORIGINAL ISSUE
DATE: 97/04/28
BY: JMH

05
1 2

2 ± 1
JMH 97/04/28
MC 97/01/5
MH 97/01/5

RF TEST CABLE

S92902X2.dgn 92902-0002 50-92902-003 1 OF 1

| ITEM | DESCRIPTION | PART NO. |
|------|---------------------------------|------------|
| ▽ | MOLEX RIGHT ANGLED RF PLUG | 92905-0001 |
| ▽ | STRAIGHT BNC CABLE (RG178) JACK | - |
| ▽ | RG178 CABLE | - |

MODIFICATION RECORD INDEX

| ECN | ISS. | DATE | INIT. | SHTS. AFFECTED | DESCRIPTION OF MOD. |
|--------|------|----------|-------|----------------|--|
| E70683 | A | 97.05.20 | JMH | 1 | ORIGINAL ISSUE |
| | | | | | |
| | | | | | REVISE ONLY ON CAD SYSTEM CHECK FILENAME: S92902X2 .DGN ISSUE: A |

LEVEL CHART

| | | | |
|----|----|----|-------------------------|
| 1 | 17 | 33 | 49 |
| 2 | 18 | 34 | 50 |
| 3 | 19 | 35 | 51 LEVEL INDEX |
| 4 | 20 | 36 | 52 RESERVED |
| 5 | 21 | 37 | 53 RESERVED |
| 6 | 22 | 38 | 54 Inspection Balloons |
| 7 | 23 | 39 | 55 DIMENSIONS |
| 8 | 24 | 40 | 56 TABLES/CHARTS |
| 9 | 25 | 41 | 57 REVISION TEXT & MRI |
| 10 | 26 | 42 | 58 Notes/text - English |
| 11 | 27 | 43 | 59 Bill of materials |
| 12 | 28 | 44 | 60 notes/text - german |
| 13 | 29 | 45 | 61 |
| 14 | 30 | 46 | 62 notes/text - french |
| 15 | 31 | 47 | 63 RESERVED |
| 16 | 32 | 48 | |