



Directions:

Release fluid from swab handle into tip by carefully bending handle. Safety gloves recommended. Discard following use.

For industrial use only. Read SDS carefully prior to use.



8020 Staticide[®] IPA Clean Swab

IPA-ready for critical cleaning

IPA Clean Swab is a high quality, cleanroom swab for the critical cleaning of small and hard-to-reach areas. The wide paddle head offers extraordinary cleaning abilities by featuring an inner layer of foam for absorbency and an outer layer of micro-denier for holding contaminants.

A four-inch handle containing a reservoir of 99.9% isopropyl alcohol adds even more cleaning power to the paddle head. Once the solvent is released from the handle, it is absorbed into the soft paddle head. This easy-to-use system is ideal for removing lubricants, adhesives, residues, or other contaminants on circuit boards, sensors, or stencils.

Class 100 laundered and packaged, IPA Clean Swab is free of particulates making it suitable for medical surfaces. With no tribogeneration, the swab is safe to use on ESD static-sensitive devices when the head is moist.

Ideal for the following applications:

- Medical device manufacturing and assembly
- Electronics manufacturing and assembly
- IoT and Telecom applications

FEATURES:

- **®** ISO Class 4-7 compatible
- Ultra-low particle fiber, NVR, and Ion levels
- **# High solvent capacity**
- Micro-denier polyester over foam
- Non-abrasive wide paddle head
- 99.7% IPA-filled 4-inch polypropylene handle
- **ℜ** Free of silicone, amides, DOP contamination
- Safe for ESD static-sensitive surfaces



Product# 8020

50 individually wrapped swabs per box; 5 boxes per case

ACL Staticide

840 W. 49th Place · Chicago, IL 60609 800.782.8420 · 847.981.9212 · fax 847.981.9278 info@aclstaticide.com · www.aclstaticide.com



ACL Inc. registered to ISO-9001 : 2015 Certificate NO. 10002746

Rev. 04/43/2019 MKB